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Intended Use

The COVID-19 direct RT-PCR **Omicron** is an *in vitro* real-time test based on the reverse transcriptase polymerase chain reaction (RT-qPCR) for the qualitative detection of nucleic acid from SARS-CoV-2 in samples of the upper and lower respiratory tract collected from suspected COVID-19 individuals by their healthcare provider. Positive results permit the **identification of SARS-CoV-2 RNA and the omicron variant, respectively (via the mutation Y505H)**. Clinical correlation with medical records and other diagnostic information is necessary to determine the patient's infection status.

Package Content

Vials with solution **A_{Omi}** (1.6 mL) and solution **B_{Omi}** (18 µL), each for one microtiter plate (96 well; not included), positive control PC (13 µL, optional), negative control NC (13 µL, optional), instructions for use.

Notes before starting

The starting material for the COVID-19 direct RT-PCR **Omicron** test is **10µL/ reaction solubilised patient sample** obtained from sputum or smear. Heat inactivation of the patient sample beforehand does not affect the performance of the reagents used in the test. A possible influence of the heat inactivation on the sample can not be excluded. One positive and one negative control should be included in each PCR run. A reaction vessel corresponds to a single Eppendorf tube, one tube of a PCR-stripe or one well of a 96-well microtiter plate.

Material provided by user

- Adequate pipettes and sterile filter-tips for PCR testing
- PCR-reaction vessel and sealing foil or lid
- qPCR instrument (e.g. LightCycler 480)

Test procedure

1. **Thaw all reagents** completely and **keep them cool (+2 °C to +8 °C)** directly before starting the test.
2. Prepare the reaction solution¹: **add 15 µL of solution B_{Omi} to the vial of solution A_{Omi}**; mix/shake briefly and centrifuge if necessary. **Do not vortex!**
3. Pipette **15 µL/reaction vessel** of the **reaction solution**.
4. **Add 10 µL/reaction vessel** of a solubilised **patient sample** or Positive Control or Negative Control.
5. Close the reaction vessel with an adhesive optical film or with the lids provided.
6. Briefly centrifuge the reaction vessel, if necessary.
7. Place the filled reaction vessel in the qPCR cycler. Start program.

Instrument settings

Steps	Temperature [°C]	Time	Number of cycles
Reverse Transcription	55	10 min	1x
Denaturation	95	2 min	1x
Denaturation	95	5 sec	45x
Amplification/Elongation ²	61	15 sec	

Detection settings

	SARS-CoV-2 (N gene)	Internal Control	SARS-general (E gene)	SARS-CoV-2 omicron mutation Y505H (spike gene)
Reporter dye	FAM probe	HEX probe	CalFluor Red 610 probe	Cy5 probe
Emission	510 nm	580 nm	610 nm	664 nm
Quencher	Black Hole Quencher			

¹Use reaction solution within 4 h after preparation. Ensure continuous storage of the reaction solution at +2 °C to +8 °C. ²Activate the detection mode of your qPCR-instrument.

Interpretation of COVID-19 direct RT-PCR Omicron test results

SARS-general (E gene)	SARS-CoV-2 (N gene)	Internal Control	Result	Interpretation
Positive	Positive	Positive/Negative	Positive!	SARS corona viruses detected in the examined material; confirmation for SARS-CoV-2. If the signal in the Cy5 channel is positive, the omicron variant (Y505H mutation) of SARS-CoV-2 is present.
Positive	Negative	Positive/Negative	Potentially positive	SARS corona virus detected in the examined material; no confirmation for SARS-CoV-2; re-analysis of sample if necessary. If the signal in the Cy5 channel is positive, the omicron variant (Y505H mutation) of SARS-CoV-2 is present.
Negative	Negative	Positive	Negative!	No SARS corona virus and no SARS-CoV-2 were found in the examined material.
Negative	Positive	Positive/Negative	Potentially positive	No SARS corona virus detected in the material examined, but SARS-CoV-2; re-analysis of sample if necessary. If the signal in the Cy5 channel is positive, the omicron variant (Y505H mutation) of SARS-CoV-2 is present.
Negative	Negative	Negative	Invalid	The test result cannot be evaluated.

The results are used to identify SARS-CoV-2 RNA. The SARS-CoV-2-RNA is generally detectable in respiratory samples during the acute phase of the infection. Positive results are an indication of the presence of SARS-CoV-2 RNA. A negative result does not rule out the presence of SARS-CoV-2 RNA, as the results depend on correct sampling and a sufficient amount of RNA to be detected.

Important notes

All samples of biological origin and used plates/swabs are to be treated as potential carriers of infectious diseases.

When working with chemicals or when handling samples of biological origin, the safety precautions of the laboratory must be observed.

Before performing this test, read the instructions for use to familiarise yourself with the testing procedure. You can find them on <https://frizbiochem.de/downloads/>.

If you have any questions or problems, please contact service at FRIZ Biochem GmbH (<https://frizbiochem.de/get-in-touch/>).

Storage: -25°C to -18°C

Usage: +2°C to +8°C

Disposal: According to the local, regional or governmental regulations

Quality control: Internal Control in the test; 1 positive control per run; 1 negative control per run

Cross reactions: None found